


Jasco V770 UV/VIS Spectrophotometer user manual

- Switch on spectrophotometer
- Boot up PC
- Start Program „Spectra Manager 2.5“

Program manual

- Double-click on „Spectra Measurement“ to open system control
- If necessary take out previous samples, if the sample is unknown contact a supervisor and do not touch until danger is excluded
Please leave minimize the opening time of the cover to prevent sensor damage
- When using a reference, insert it nor, if necessary in both spots
- Press „Baseline“
- Insert sample einstecken (sample in front, reference in back)
- Press „Sample“

*Change parameter: „Measure“
„Parameter“ and adjust by your liking:*



Examine and extract data

- You can view your data in the window „Spectra Analysis“ and save it (recommended is “Collectively CSV”) after having chosen the wanted curve “File” -> “Export”
The instrument can't measure an absorption of 5 or higher, it will jump to 10, don't get confused

End of experiment

- End program
- Switch off spectrometer
- Take out samples & references
- Close the lid